


Anti-Pghs-2 Polyclonal Antibody For Wb - P35354

Item Number: CM0019801



Introduction

Rabbit polyclonal antibody targeting human PGHS-2 (COX-2/PTGS2), suitable for Western blot. This enzyme plays a central role in prostanoid biosynthesis and inflammatory response. Use for inflammation and lipid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	PGHS-2 Rabbit pAb
Remarks/Alias	COX2; COX-2; PHS-2; PGG/HS; PGHS-2; hCox-2; GRIPGHS
Species Reactivity	Human
GeneID (Human)	5743
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 505-604 of human COX2/PTGS2 (P35354).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB
Cross Reactivity	Human
SWISS	P35354
Protein Weight	69kDa
Shipping	Ice bag