

Anti-Pgk1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P00558

Item Number: CM0014254



Introduction

Rabbit polyclonal antibody against PGK1 (Phosphoglycerate kinase 1). Validated in WB, IHC-P, ELISA; cross-reacts with human, mouse, rat. For research on glycolysis, mitochondrial function, and cell migration.

[Learn More](#)

Parameter	Value
Product Name	PGK1 Rabbit pAb
Remarks/Alias	PGKA; MIG10; HEL-S-68p; PGK1
Species	Human
GeneID (Human)	5230
GeneID	5230
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-417 of human PGK1 (NP_000282.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00558
Protein Weight	45kDa
Shipping	Ice bag