

Anti-Pgk2 Rabbit Polyclonal Antibody For Wb, Ihc, If/Icc, Elisa - P07205

Item Number: CM0015924

Introduction

Rabbit polyclonal PGK2 antibody validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat PGK2 (45 kDa), a testis-specific glycolytic enzyme essential for sperm motility.

[Learn More](#)

Parameter	Value
Product Name	PGK2 Rabbit pAb
Remarks/Alias	PGKB; PGKPS; HEL-S-272; dj417L20.2; PGK2
Species	Human
GeneID (Human)	5232
GeneID	5232
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 242-360 of human PGK2 (NP_620061.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07205
Protein Weight	45kDa
Shipping	Ice bag