

Anti-Pgls Polyclonal Antibody For Wb, Ihc-P, Elisa - 095336

Item Number: CM0028701

Introduction

Rabbit polyclonal antibody targeting human 6-phosphogluconolactonase (PGLS). Validated for WB, IHC-P, ELISA with cross-reactivity to mouse and rat. Ideal for metabolic pathway studies.

[Learn More](#)

Parameter	Value
Product Name	PGLS Rabbit pAb
Remarks/Alias	6PGL; HEL-S-304; PGLS
Species	Human
GeneID	25796
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-258 of human PGLS (NP_036220.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95336
Protein Weight	28kDa
Shipping	Ice bag