

Anti-Pglyrp2 Polyclonal Antibody For Wb, Elisa - Q96Pd5

Item Number: CM0021071

Introduction

Rabbit polyclonal antibody against human PGLYRP2 (peptidoglycan recognition protein 2, PGRP-L), suitable for WB and ELISA. Cross-reacts with mouse. Ideal for innate immunity and peptidoglycan digestion research.

[Learn More](#)

Parameter	Value
Product Name	PGLYRP2 Rabbit pAb
Remarks/Alias	tagL; PGRPL; PGRP-L; PGLYRPL; HMFT0141; TAGL-like; tagl-beta; tagl-alpha; PGLYRP2
Species	Human
Gene ID (Human)	114770
Gene ID	114770
Immunogen	Synthetic peptide (A synthetic peptide corresponding to a sequence within amino acids 300-400 of human PGLYRP2 (NP_443122.3))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q96PD5
Protein Weight	62kDa
Shipping	Ice bag