

Anti-Pgm2 Rabbit Polyclonal Antibody For Wb, Elisa - Q96G03

Item Number: CM0015217

Introduction

PGM2 Rabbit Polyclonal Antibody is a polyclonal antibody raised in rabbit against human PGM2, validated for WB and ELISA. Detects human and mouse PGM2, a phosphopentomutase crucial for nucleotide metabolism. Shipped on ice.

[Learn More](#)

Parameter	Specification
Product Name	PGM2 Rabbit pAb
Remarks/Alias	MSTP006; PGM2
Species	Human
Gene ID (Human)	55276
Gene ID	55276
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human PGM2 (NP_060760.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96G03
Protein Weight	68kDa
Shipping	Ice bag