


Anti-Pgp9.5/Uchl1 Polyclonal Antibody For Wb, Elisa - P09936

Item Number: CM0018096



Introduction

High-quality rabbit polyclonal antibody targeting PGP9.5/Uchl1 (UniProt P09936). Suitable for Western blot (WB) and ELISA applications. Cross-reacts with human and mouse Uchl1, a key deubiquitinase implicated in neurodegeneration and protein homeostasis.

[Learn More](#)

Parameter	Value
Product Name	PGP9.5/Uchl1 Rabbit pAb
Remarks/Alias	NDGOA; PARK5; PGP95; SPG79; PGP9.5; SPG79A; Uchl-1; Uch-L1; HEL-117; PGP 9.5; HEL-S-53; PGP9.5/Uchl1
Species	Human
GeneID (Human)	7345
GeneID	7345
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P09936
Protein Weight	25kDa
Shipping	Ice bag