


Anti-Pgp9.5/Uchl1 Polyclonal Antibody For Wb, If/Icc, Elisa - P09936

Item Number: CM0021729



Introduction

Rabbit polyclonal antibody against human PGP9.5/UCHL1 (ubiquitin carboxyl-terminal hydrolase isozyme L1), a key deubiquitinase in protein degradation and neuronal function. Validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for neuroscience and neurodegenerative disease research.

[Learn More](#)

Parameter	Value
Product Name	PGP9.5/UCHL1 Rabbit pAb
Remarks/Alias	NDGOA; PARK5; PGP95; SPG79; PGP9.5; SPG79A; UCHL-1; Uch-L1; HEL-117; PGP 9.5; HEL-S-53; PGP9.5/UCHL1
Species	Human
GeneID (Human)	7345
GeneID	7345
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human PGP95/UCHL1 (NP_004172.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09936
Protein Weight	25kDa
Shipping	Ice bag