


Anti-Pgp9.5/Uchl1 Rabbit Monoclonal Antibody For Wb, If, Ihc, Elisa - P09936

Item Number: CM0002236



Introduction

Rabbit monoclonal antibody to human PGP9.5/UCHL1 (25 kDa), validated for WB, IF, IHC, ELISA; reacts with human, mouse, rat. A deubiquitinase for synaptic and cardiac function, APP processing, osteoclastogenesis, and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	PGP9.5/UCHL1 Rabbit mAb
Remarks/Alias	NDGOA; PARK5; PGP95; SPG79; PGP9.5; SPG79A; UCHL-1; Uch-L1; HEL-117; PGP 9.5; HEL-S-53; PGP9.5/UCHL1
Species	Human
GeneID (Human)	7345
GeneID	7345
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human PGP95/UCHL1 (NP_004172.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09936
Protein Weight	25kDa
Shipping	Ice bag