


## Anti-Pgp9.5/Uchl1 Rabbit Monoclonal Antibody For Wb, If/Icc, If-F, If-P, Ihc-P, Elisa - P09936

Item Number: CM0003461



### Introduction

Rabbit monoclonal antibody against human PGP9.5/UCHL1 (SWISS P09936). Validated for WB, IF/ICC, IHC-P, ELISA; cross-reacts with human, mouse, rat, monkey. Suitable for neurodegeneration and ubiquitin pathway research.

[Learn More](#)

| Parameter        | Value                                                                                                                    |
|------------------|--------------------------------------------------------------------------------------------------------------------------|
| Product Name     | PGP9.5/UCHL1 Rabbit mAb                                                                                                  |
| Remarks/Alias    | NDGOA; PARK5; PGP95; SPG79; PGP9.5; SPG79A; UCHL-1; Uch-L1; HEL-117; PGP 9.5; HEL-S-53; PGP9.5/UCHL1                     |
| Species          | Human                                                                                                                    |
| GeneID (Human)   | 7345                                                                                                                     |
| GeneID           | 7345                                                                                                                     |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human PGP9.5/UCHL1 (NP_004172.2). |
| Source           | Rabbit                                                                                                                   |
| Category         | Monoclonal Antibodies                                                                                                    |
| Application      | WB, IF/ICC, IF-F, IF-P, IHC-P, ELISA                                                                                     |
| Cross Reactivity | Human, Mouse, Rat, Monkey                                                                                                |
| SWISS            | P09936                                                                                                                   |
| Protein Weight   | 25 kDa                                                                                                                   |
| Shipping         | Ice bag                                                                                                                  |