


Anti-Pgrmc1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O00264

Item Number: CM0028362



Introduction

Rabbit polyclonal antibody against human PGRMC1, validated for WB, IHC-P, IF/ICC, and ELISA. Detects endogenous levels in human, mouse, and rat. Suitable for progesterone signaling and heme homeostasis research.

[Learn More](#)

Parameter	Value
Product Name	PGRMC1 Rabbit pAb
Remarks/Alias	IZA; MPR; Dap1; HPR6.6; PGRMC1
Species	Human
GeneID (Human)	10857
GeneID	10857
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 44-195 of human PGRMC1 (NP_006658.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00264
Protein Weight	22kDa
Shipping	Ice bag