


Anti-Phactr1 Polyclonal Antibody For Wb, Elisa - Q9C0D0

Item Number: CM0023140



Introduction

CamelBio offers Anti-PHACTR1 Rabbit Polyclonal Antibody for WB and ELISA. Recognizes human, mouse, and rat PHACTR1, a phosphatase and actin regulator implicated in cell motility and cytoskeletal organization.

[Learn More](#)

Parameter	Value
Product Name	PHACTR1 Rabbit pAb
Remarks/Alias	RPEL; DEE70; RPEL1; EIEE70; dJ257A7.2; PHACTR1
Species	Human
GeneID (human)	221692
GeneID	221692
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 154-380 of human PHACTR1 (NP_112210.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9C0D0
Protein Weight	66kDa
Shipping	Ice bag