

Anti-Phc3 Polyclonal Antibody For Wb, Elisa - Q8Ndx5

Item Number: CM0028748



Introduction

Rabbit polyclonal antibody targeting human PHC3, validated for Western blot and ELISA. Reacts with human, mouse, and rat. Suitable for Polycomb group protein and chromatin remodeling research.

[Learn More](#)

Parameter	Value
Product Name	PHC3 Rabbit pAb
Remarks/Alias	EDR3; HPH3; PHC3
Species	Human
GeneID (Human)	80012
GeneID	80012
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-300 of human PHC3 (NP_079223.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NDX5
Protein Weight	106kDa
Shipping	Ice bag