


Anti-Phc3 Rabbit Polyclonal Antibody For If/Icc, Elisa - Q8Ndx5

Item Number: CM0017458



Introduction

Rabbit polyclonal antibody against human PHC3, validated for IF/ICC and ELISA. Targets the polyhomeotic-like protein 3, a component of the PRC1 complex, for epigenetic and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	PHC3 Rabbit pAb
Remarks/Alias	EDR3, HPH3, PHC3
Species	Human
GeneID (Human)	80012
GeneID	80012
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-300 of human PHC3 (NP_079223.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q8NDX5
Protein Weight	106kDa
Shipping	Ice bag