

Anti-Phd2/Egln1 Polyclonal Antibody For If/Icc, Elisa - Q9Gzt9

Item Number: CM0024305



Introduction

Rabbit polyclonal antibody targeting PHD2/EGLN1 (HIF prolyl hydroxylase 2), validated for IF/ICC and ELISA. Detects human PHD2, a key oxygen sensor regulating HIF-alpha stability. Ideal for hypoxia and angiogenesis research.

[Learn More](#)

Parameter	Value
Product Name	PHD2/EGLN1 Rabbit pAb
Remarks/Alias	HPH2; PHD2; SM20; ECTY3; HALAH; HPH-2; HIFPH2; ZMYND6; C1orf12; HIF-PH2; PHD2/EGLN1
Species	Human
GeneID (human)	54583
GeneID	54583
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-426 of human PHD2/EGLN1 (NP_071334.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9GZT9
Protein Weight	46kDa
Shipping	Ice bag