

Anti-Phd2/Egln1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Gzt9

Item Number: CM0018045



Introduction

Rabbit polyclonal antibody against human PHD2/EGLN1, validated for WB, IHC-P, IF/ICC, and ELISA. Detects endogenous PHD2/EGLN1 in human, mouse, and rat samples. Ideal for hypoxia and angiogenesis research.

[Learn More](#)

Parameter	Value
Product Name	PHD2/EGLN1 Rabbit pAb
Remarks/Alias	HPH2; PHD2; SM20; ECTY3; HALAH; HPH-2; HIFPH2; ZMYND6; C1orf12; HIF-PH2; PHD2/EGLN1
Species	Human
GeneID (Human)	54583
GeneID	54583
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-426 of human PHD2/EGLN1 (NP_071334.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZT9
Protein Weight	46kDa
Shipping	Ice bag