

Anti-Phf1 Rabbit Polyclonal Antibody For Wb, Elisa - O43189

Item Number: CM0027630



Introduction

High-quality anti-PHF1 rabbit polyclonal antibody validated for Western blot and ELISA. Detects human, mouse, and rat PHF1. Ideal for chromatin research, DNA damage response, and epigenetic studies.

[Learn More](#)

Parameter	Value
Product Name	PHF1 Rabbit pAb
Remarks/Alias	PCL1; PHF2; hPHF1; MTF2L2; TDRD19C; PHF1
Species	Human
GeneID (Human)	5252
GeneID	5252
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human PHF1 (NP_077084.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43189
Protein Weight	62kDa
Shipping	Ice bag