

Anti-Phf19 Polyclonal Antibody For Wb, If/Icc, Elisa - Q5T6S3

Item Number: CM0028611



Introduction

Rabbit polyclonal antibody targeting PHF19/PCL3, validated for WB, IF/ICC, ELISA applications. Detects human, mouse, and rat PHF19. Ideal for studying Polycomb repression, H3K36me3 recognition, and cancer epigenetics.

[Learn More](#)

Parameter	Value
Product Name	PHF19 Rabbit pAb
Remarks/Alias	PCL3; MTF2L1; TDRD19B; PHF19
Species	Human
GeneID (Human)	26147
GeneID	26147
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-207 of human PHF19 (NP_001009936.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5T6S3
Protein Weight	66kDa
Shipping	Ice bag