

Anti-Phf21B Rabbit Polyclonal Antibody For Wb & Elisa - Q96Ek2

Item Number: CM0028472

Introduction

Rabbit polyclonal antibody targeting PHF21B (PHF4, BHC80L), validated for WB and ELISA in human, mouse, and rat samples. Ideal for chromatin biology and epigenetics research.

[Learn More](#)

Parameter	Value
Product Name	PHF21B Rabbit pAb
Remarks/Alias	PHF4; BHC80L; PHF21B
Species	Human
GeneID (Human)	112885
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 322-531 of human PHF21B (NP_612424.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96EK2
Protein Weight	57kDa
Shipping	Ice bag