

Anti-Phf8 Monoclonal Antibody For Wb, Elisa - Q9Upp1

Item Number: CM0002912



Introduction

Rabbit monoclonal antibody against PHF8 (KDM7B), histone lysine demethylase, validated for WB and ELISA in human and mouse samples. Ideal for epigenetic and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	PHF8 Rabbit mAb
Remarks/Alias	KDM7B; JHDM1F; MRXSSD; ZNF422; PHF8
Species	Human
GeneID (human)	23133
GeneID	23133
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 800-900 of human PHF8 (Q9UPP1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UPP1
Protein Weight	118kDa
Shipping	Ice bag