

Anti-Phgdh Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O43175

Item Number: CM0004220



Introduction

Recombinant rabbit monoclonal antibody against human PHGDH (O43175). Validated for WB, IHC-P, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for studying serine biosynthesis and cancer metabolism.

[Learn More](#)

Parameter	Value
Product Name	PHGDH Rabbit mAb
Remarks/Alias	NLS; PDG; PGD; NLS1; PGAD; PGDH; SERA; 3PGDH; 3-PGDH; PHGDHD; HEL-S-113; PHGDH
Species	Human
Gene ID (Human)	26227
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 264-533 of human PHGDH (NP_006614.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43175
Protein Weight	57kDa
Shipping	Ice bag