

Anti-Phgdh Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - O43175

Item Number: CM0017947



Introduction

Rabbit polyclonal antibody recognizing PHGDH, a key enzyme in serine biosynthesis. Validated for WB, IF/ICC, IP, ELISA in human, mouse, rat. Ideal for diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	PHGDH Rabbit pAb
Remarks/Alias	NLS; PDG; PGD; NLS1; PGAD; PGDH; SERA; 3PGDH; 3-PGDH; PHGDHD; HEL-S-113; PHGDH
Species	Human
GeneID (Human)	26227
GeneID	26227
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 264-533 of human PHGDH (NP_006614.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43175
Protein Weight	57kDa
Shipping	Ice bag