


Anti-Phip Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q8Wwq0

Item Number: CM0024788



Introduction

Rabbit polyclonal antibody targeting human PHIP (PH-interacting protein), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse. Ideal for insulin signaling and cell proliferation studies.

[Learn More](#)

Parameter	Value
Product Name	PHIP Rabbit pAb
Remarks/Alias	ndrp; BRWD2; DIDOD; WDR11; DCAF14; CHUJANS; PHIP
Species	Human
GeneID (Human)	55023
GeneID	55023
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1542-1821 of human PHIP (NP_060404.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8WWQ0
Protein Weight	207kDa
Shipping	Ice bag