


Anti-Phka1 Polyclonal Antibody For Wb, Elisa - P46020

Item Number: CM0021552



Introduction

PHKA1 Rabbit Polyclonal Antibody validated for WB and ELISA. Recognizes human and mouse phosphorylase b kinase regulatory subunit alpha, ideal for glycogen metabolism and muscle disease research.

[Learn More](#)

Parameter	Value
Product Name	PHKA1 Rabbit pAb
Remarks/Alias	PHKA; PHKA1
Species	Human
GeneID	5255
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 600-820 of human PHKA1 (NP_001116142.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P46020
Protein Weight	137kDa
Shipping	Ice bag