


Anti-Phkb Polyclonal Antibody For Wb, Elisa - Q93100

Item Number: CM0026054



Introduction

Rabbit polyclonal antibody targeting human PHKB, the regulatory subunit beta of phosphorylase b kinase. Suitable for WB and ELISA, with cross-reactivity to human and mouse. Ideal for glycogen metabolism and related disease studies.

[Learn More](#)

Parameter	Value
Product Name	PHKB Rabbit pAb
Remarks/Alias	PHKB
Species	Human
Gene ID (Human)	5257
Gene ID	5257
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 854-1093 of human PHKB (NP_000284.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q93100
Protein Weight	125kDa
Shipping	Ice bag