


Anti-Phkg2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P15735

Item Number: CM0025716



Introduction

Rabbit polyclonal antibody targeting human PHKG2 (phosphorylase b kinase gamma 2). Suitable for WB, IF/ICC, and ELISA with cross-reactivity to mouse and rat. Research use in glycogenolysis and GSD type IX studies.

[Learn More](#)

Parameter	Value
Product Name	PHKG2 Rabbit pAb
Remarks/Alias	GSD9C; PHKG2
Species	Human
GeneID (Human)	5261
GeneID	5261
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 237-406 of human PHKG2 (NP_000285.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15735
Protein Weight	46kDa
Shipping	Ice bag