

Anti-Phlpp1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O60346

Item Number: CM0025054

Introduction

Rabbit polyclonal antibody targeting PHLPP1, validated for WB, IHC-P, IF/ICC, ELISA. Recognizes human, mouse, and rat PHLPP1. Ideal for research on apoptosis, tumor suppression, and circadian biology.

[Learn More](#)

Parameter	Value
Product Name	PHLPP1 Rabbit pAb
Remarks/Alias	SCOP; PHLPP; PPM3A; PLEKHE1; PHLPP1
Species	Human
GeneID (Human)	23239
GeneID	23239
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1160-1428 of human PHLPP1 (NP_919431.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60346
Protein Weight	185kDa
Shipping	Ice bag