

Anti-Phlpp2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q6Zvd8

Item Number: CM0020597

Introduction

Rabbit polyclonal antibody against human PHLPP2 (PHLPP-like), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Target is a tumor suppressor phosphatase regulating Akt and PKC signaling, with roles in apoptosis and proliferation.

[Learn More](#)

Parameter	Value
Product Name	PHLPP2 Rabbit pAb
Remarks/Alias	PPM3B; PHLPL; PHLPP2
Species	Human
Gene ID (Human)	23035
Gene ID	23035
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1250-1323 of human PHLPP2 (NP_055835.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6ZVD8
Protein Weight	147kDa
Shipping	Ice bag