

Anti-Phospho-Akt-T308 Rabbit Monoclonal Antibody For Wb, Elisa - P31749 / P31751 / Q9Y243

Item Number: CM0000931



Introduction

High-specificity phospho-Akt-T308 rabbit monoclonal antibody for WB and ELISA. Recognizes human Akt1/2/3 phosphorylated at T308; cross-reacts with mouse and rat. Ideal for studying PI3K/Akt signaling in cancer and metabolism research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Akt-T308 Rabbit mAb
Remarks/Alias	AKT1/AKT2/AKT3; Phospho-Akt-T308
Species	Human
Gene ID (Human)	207/208/10000
Gene ID	207/208/10000
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around T308 of human Akt (NP_005154.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P31749/P31751/Q9Y243
Protein Weight	48kDa/55kDa/51kDa/54kDa
Shipping	Ice bag