

Anti-Phospho-Akt1-S473 Monoclonal Antibody For Wb, Ihc-P, Elisa - P31749

Item Number: CM0001091



Introduction

Rabbit monoclonal antibody against phosphorylated Akt1 at Ser473 (P31749). Validated for WB, IHC-P, ELISA in human, mouse, rat samples. Ideal for studying Akt signaling in cancer and metabolism.

[Learn More](#)

Parameter	Value
Product Name	Phospho-AKT-S473 Rabbit mAb
Remarks/Alias	AKT; PKB; RAC; PRKBA; PKB-ALPHA; RAC-ALPHA; Phospho-AKT1-S473
Species	Human
GeneID (Human)	207
GeneID	207
Immunogen	A synthetic phosphorylated peptide around S473 of human Akt1 (P31749)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31749
Protein Weight	56kDa
Shipping	Ice bag