

Anti-Phospho-Akt1 (Ser473) Rabbit Polyclonal Antibody For Wb, Elisa - P31749

Item Number: CM0001196



Introduction

Phospho-AKT1 (Ser473) rabbit polyclonal antibody validated for WB and ELISA. Detects human, mouse, and rat AKT1 phosphorylated at Ser473. Suitable for studying AKT signaling in cancer, metabolism, and cell survival. Sourced from immunized rabbits.

[Learn More](#)

Parameter	Value
Product Name	Phospho-AKT1-S473 Rabbit pAb
Remarks/Alias	AKT; PKB; RAC; PRKBA; PKB-ALPHA; RAC-ALPHA; Phospho-AKT1-S473
Species	Human
GeneID (Human)	207
GeneID	207
Immunogen	Synthetic peptide; A synthetic phosphorylated peptide around S473 of human Akt1 (NP_005154.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31749
Protein Weight	56kDa
Shipping	Ice bag