

Anti-Phospho-Akt1-T308 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P31749

Item Number: CM0001741



Introduction

Phospho-specific rabbit polyclonal antibody targeting AKT1 phosphorylated at Thr-308. Validated for Western Blot, IHC-P, and ELISA. Recognizes human, mouse, and rat AKT1. Ideal for studying PI3K/AKT signaling in cancer and metabolism research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-AKT1-T308 Rabbit pAb
Remarks/Alias	AKT; PKB; RAC; PRKBA; PKB-ALPHA; RAC-ALPHA; Phospho-AKT1-T308
Species	Human
GeneID (Human)	207
GeneID	207
Immunogen	Synthetic peptide – A synthetic phosphorylated peptide around T308 of human AKT (NP_005154.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31749
Protein Weight	56kDa
Shipping	Ice bag