

Anti-Phospho-Akt1 (T450) Rabbit Monoclonal Antibody For Wb, Elisa - P31749

Item Number: CM0000681



Introduction

Phospho-AKT1-T450 Rabbit mAb for WB and ELISA; specifically detects AKT1 phosphorylated at Thr450 in human and rat samples. Ideal for studying PI3K/AKT signaling in cancer, metabolism, and cell survival.

[Learn More](#)

Parameter	Value
Product Name	Phospho-AKT1-T450 Rabbit mAb
Remarks/Alias	AKT; PKB; RAC; PRKBA; PKB-ALPHA; RAC-ALPHA; Phospho-AKT1-T450
Species	Human
GeneID (Human)	207
GeneID	207
Immunogen	A synthetic phosphorylated peptide around T450 of human Akt1 (P31749).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P31749
Protein Weight	56kDa
Shipping	Ice bag