


Anti-Phospho-Alox5-S271 Polyclonal Antibody For Wb, Elisa - P09917

Item Number: CM0001771



Introduction

Rabbit polyclonal antibody specific to human ALOX5 phosphorylated at serine 271, validated for Western blot and ELISA applications in leukotriene biosynthesis and inflammatory disease research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ALOX5-S271 Rabbit pAb
Remarks/Alias	5-LO; 5LPG; LOG5; 5-LOX; Phospho-ALOX5-S271
Species	Human
Gene ID (Human)	240
Gene ID	240
Immunogen	A synthetic phosphorylated peptide around S271 of human ALOX5 (NP_000689.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P09917
Protein Weight	78kDa
Shipping	Ice bag