


Anti-Phospho-Atf2 (T69/T71) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P15336

Item Number: CM0000897



Introduction

Rabbit polyclonal antibody targeting phosphorylated human ATF2 at Thr69/Thr71. Validated for WB, IHC-P, and ELISA; cross-reacts with mouse and rat. Ideal for DNA damage response, apoptosis, and transcription studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ATF2-T69/71 Rabbit pAb
Remarks/Alias	HB16; CREB2; TREB7; CREB-2; CRE-BP1; Phospho-ATF2-T69/71
Species	Human
Gene ID (Human)	1386
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15336
Protein Weight	55kDa
Shipping	Ice bag