

Anti-Phospho-Atm (S1981) Rabbit Polyclonal Antibody For Wb And Elisa - Q13315

Item Number: CM0001396



Introduction

Rabbit polyclonal antibody specific for ATM phosphorylated at Ser1981. Validated for WB and ELISA in human and mouse samples. Ideal for DNA damage signaling and cell cycle checkpoint studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ATM-S1981 Rabbit pAb
Remarks/Alias	AT1; ATA; ATC; ATD; ATE; ATDC; TEL1; TEL01; Phospho-ATM-S1981
Species	Human
Gene ID (Human)	472
Gene ID	472
Immunogen	A synthetic phosphorylated peptide around S1981 of human ATM (NP_000042.3)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13315
Protein Weight	351 kDa
Shipping	Ice bag