


Anti-Phospho-Axl (Y702) Rabbit Monoclonal Antibody For Wb, Elisa - P30530

Item Number: CM0001926



Introduction

Phospho-specific rabbit monoclonal antibody targeting human Axl phosphorylated at Tyr702. Validated for WB and ELISA. Ideal for studying Axl kinase signaling, cell survival, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Axl-Y702 Rabbit mAb
Remarks/Alias	ARK; UFO; JTK11; Tyro7; Phospho-Axl-Y702
Species	Human
GeneID	558 (Human)
Immunogen	A synthetic phosphorylated peptide around Y702 of human Axl (NP_068713.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P30530
Protein Weight	98 kDa
Shipping	Ice bag