

Anti-Phospho-B-Raf (T401) Rabbit Monoclonal Antibody For Wb, Elisa - P15056

Item Number: CM0001347



Introduction

High-quality rabbit monoclonal antibody targeting phosphorylated B-Raf at Thr401. Suitable for WB and ELISA with rat cross-reactivity. Investigate MAP kinase signaling and oncogenic pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-B-Raf-T401 Rabbit mAb
Remarks/Alias	NS7; B-raf; BRAF1; RAFB1; B-RAF1; BRAF-1; Phospho-B-Raf-T401
Species	Human
GeneID (Human)	673
GeneID	673
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P15056
Protein Weight	84kDa
Shipping	Ice bag