


Anti-Phospho-Bad-S112 Rabbit Monoclonal Antibody For Wb And Elisa - Q92934

Item Number: CM0000978



Introduction

Phospho-Bad-S112 rabbit monoclonal antibody validated for WB and ELISA, detecting human and mouse Bad phosphorylated at serine 112. 18 kDa. Ideal for apoptosis research. From CamelBio, a one-stop IVD raw material provider.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Bad-S112 Rabbit mAb
Remarks/Alias	BBC2; BCL2L8; Phospho-Bad-S112
Species	Human
Gene ID (Human)	572
Gene ID	572
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around S112 of human BAD (Q92934)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92934
Protein Weight	18kDa
Shipping	Ice bag