

Anti-Phospho-Bcar1 (Y410) Polyclonal Antibody For Wb - P56945

Item Number: CM0001193



Introduction

Rabbit polyclonal antibody specific to phosphorylated BCAR1 at tyrosine 410. Validated for WB and ELISA in human and mouse samples. Ideal for studying cell adhesion and migration signaling.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-BCAR1-Y410 Rabbit pAb
Remarks/Alias	CAS; CAS1; CASS1; CRKAS; P130Cas; Phospho-BCAR1-Y410
Species	Human
GeneID (Human)	9564
GeneID	9564
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around Y410 of human BCAR1 (NP_055382.2?).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P56945
Protein Weight	93kDa
Shipping	Ice bag