

Anti-Phospho-Bcl-2-S70 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P10415

Item Number: CM0000807



Introduction

Rabbit polyclonal antibody targeting phospho-Bcl-2 (Ser70). Validated for WB, IF/ICC, ELISA; reacts with human and rat Bcl-2. Suitable for apoptosis, autophagy, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Bcl-2-S70 Rabbit pAb
Remarks/Alias	Bcl-2; PPP1R50; Phospho-Bcl-2-S70
Species	Human
Gene ID	596
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S70 of human Bcl-2 (NP_000624.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	P10415
Protein Weight	26kDa
Shipping	Ice bag