

Anti-Phospho-Bcl-2 (Thr56) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P10415

Item Number: CM0001400



Introduction

High-quality rabbit polyclonal antibody specific for phosphorylated Bcl-2 (Thr56). Detects endogenous levels in human, mouse, and rat samples via WB, IHC-P, and ELISA. Essential for apoptosis and autophagy studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Bcl-2-T56 Rabbit pAb
Remarks/Alias	Bcl-2; PPP1R50; Phospho-Bcl-2-T56
Species	Human
Gene ID (Human)	596
Gene ID	596
Immunogen	A synthetic phosphorylated peptide around T56 of human Bcl-2 (NP_000624.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10415
Protein Weight	26kDa
Shipping	Ice bag