

Anti-Phospho-Bcl-XL (S62) Rabbit Polyclonal Antibody For Wb, If/Icc - Q07817

Item Number: CM0001524



Introduction

This rabbit polyclonal antibody detects phosphorylated Bcl-XL (S62) in human, mouse, and rat samples. Ideal for WB and IF/ICC applications in apoptosis and cell death studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Bcl-XL-S62 Rabbit pAb
Remarks/Alias	BCLX; BCL2L; Bcl-X; PPP1R52; BCL-XL/S; Phospho-Bcl-XL-S62
Species	Human
Gene ID (human)	598
Gene ID	598
Immunogen	A phospho specific peptide corresponding to residues surrounding S62 of human Bcl-XL
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07817
Protein Weight	26kDa
Shipping	Ice bag