


Anti-Phospho-Bcr (Y177) Rabbit Polyclonal Antibody For Wb, Elisa - P11274

Item Number: CM0001991



Introduction

Rabbit polyclonal antibody against human Phospho-BCR (Y177) validated for WB and ELISA, cross-reacts with mouse. Detects a band at 143 kDa. Ideal for studying BCR signaling, GTPase regulation, and synaptic function. - P11274

[Learn More](#)

Parameter	Value
Product Name	Phospho-BCR-Y177 Rabbit pAb
Remarks/Alias	ALL; CML; PHL; BCR1; D22S11; D22S662; Phospho-BCR-Y177
Species	Human
GeneID (Human)	613
GeneID	613
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around Y177 of human BCR (NP_004318.3).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P11274
Protein Weight	143kDa
Shipping	Ice bag