

Anti-Phospho-Brca1 (S988) Rabbit Polyclonal Antibody For Ihc-P, If/Icc, Elisa - P38398

Item Number: CM0000683



Introduction

Phospho-BRCA1-S988 rabbit polyclonal antibody validated for IHC-P, IF/ICC, ELISA. Detects human BRCA1 phosphorylated at serine 988, a key DNA repair and tumor suppressor protein. Ideal for cancer research and assay development.

[Learn More](#)

Parameter	Value
Product Name	Phospho-BRCA1-S988 Rabbit pAb
Remarks/Alias	IRIS; PSCP; BRCAI; BRCC1; FANCS; PNCA4; RNF53; BROVCA1; PPP1R53; Phospho-BRCA1-S988
Species	Human
GeneID (Human)	672
GeneID	672
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P38398
Protein Weight	208kDa
Shipping	Ice bag