


Anti-Phospho-Brd4-S484/S488 Rabbit Polyclonal Antibody For Wb, Db, Elisa - O60885

Item Number: CM0000593



Introduction

Phospho-specific rabbit polyclonal antibody targeting BRD4 phosphorylated at Ser484/Ser488. Suitable for WB, DB, and ELISA. Ideal for studying epigenetic memory, transcription regulation, and chromatin dynamics. Reacts with human BRD4.

[Learn More](#)

Parameter	Value
Product Name	Phospho-BRD4-S484/S488 Rabbit pAb
Remarks/Alias	CAP; MCAP; HUNK1; HUNKI; FSHRG4; Phospho-BRD4-S484/S488
Species	Human
Gene ID	23476
Immunogen	Synthetic peptide (A synthetic phosphorylated peptide around S484 & S488 of human BRD4 (NP_490597.1))
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, DB, ELISA
Cross Reactivity	Human
SWISS	O60885
Protein Weight	152kDa
Shipping	Ice bag