


Anti-Phospho-Brd4 (Ser1070) Rabbit Polyclonal Antibody For Wb, Elisa - O60885

Item Number: CM0001789



Introduction

Phospho-specific rabbit polyclonal antibody against BRD4 (pS1070) for WB and ELISA applications. Human reactive; ideal for studies on chromatin reading, transcriptional regulation, and DNA damage response.

[Learn More](#)

Parameter	Value
Product Name	Phospho-BRD4-S1070 Rabbit pAb
Remarks/Alias	CAP; MCAP; HUNK1; HUNKI; FSHRG4; Phospho-BRD4-S1070
Species	Human
Gene ID (Human)	23476
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S1070 of human BRD4 (NP_490597.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O60885
Protein Weight	152kDa
Shipping	Ice bag