


Anti-Phospho-Btk (Tyr223) Rabbit Monoclonal Antibody For Wb And Elisa - Q06187

Item Number: CM0001764



Introduction

Rabbit monoclonal antibody specific for phosphorylated BTK (Tyr223) validated for WB and ELISA. Ideal for studying B-cell receptor, TLR signaling, and BTK kinase activity in human samples.

[Learn More](#)

Parameter	Value
Product Name	Phospho-BTK-Y223 Rabbit mAb
Remarks/Alias	AT; ATK; BPK; XLA; IMD1; AGMX1; IGHD3; PSCTK1; Phospho-BTK-Y223
Species	Human
GeneID (Human)	695
GeneID	695
Immunogen	A synthetic phosphorylated peptide around Y223 of human BTK (Q06187)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q06187
Protein Weight	76kDa
Shipping	Ice bag