


Anti-Phospho-C-Abl-Y245 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P00519

Item Number: CM0000571



Introduction

Phospho-specific ABL1 (c-Abl) rabbit polyclonal antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat ABL1 at ~123 kDa. Suitable for kinase signaling, DNA damage, and leukemia research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-c-Abl-Y245 Rabbit pAb
Remarks/Alias	ABL; JTK7; p150; c-ABL; v-abl; CHDSKM; c-ABL1; BCR-ABL; bcr/abl; Phospho-c-Abl-Y245
Species	Human
GeneID	25
Immunogen	A synthetic phosphorylated peptide around Y245 of human ABL1 (NP_005148.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00519
Protein Weight	123kDa
Shipping	Ice bag