


Anti-Phospho-C-Jun-S63 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P05412

Item Number: CM0001982



Introduction

High-specificity rabbit monoclonal antibody against phosphorylated c-Jun (S63). Validated for WB, IF/ICC, ELISA; reacts with human, mouse, rat. Ideal for studying AP-1 transcription factor signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-c-Jun-S63 Rabbit mAb
Remarks/Alias	AP1; p39; AP-1; cJUN; c-Jun; Phospho-c-Jun-S63
Species	Human
Gene ID (Human)	3725
Gene ID	3725
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around S63 of human c-Jun (P05412).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05412
Protein Weight	36 kDa
Shipping	Ice bag